Rx Only

Please read this Operation Manual carefully
before use, and file for future reference.

Surgical Handpiece SGS
Surgical Handpiece SGA

SGS / SGA

SGS-ES / SGS-E2S
SGA-ES / SGA-E2S

MADE IN JAPAN OM-SH0102E 001

1. User and Indications for Use

User: Licensed dentists / oral surgeons

Indications for Use:

The Surgical Straight / Angle Handpiece Series is intended for use by licensed professionals for oral surgical and dental use.

Typical uses include sectioning and removal of wisdom teeth, cutting teeth, oral surgery (such as root tip resection, bone removal, bone shaping) and implant treatment (such as drilling the maxilla and mandible).

2. Precautions for Handling and Operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

WARNING

- Clean, disinfect, lubricate and sterilize the handpiece immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the handpiece may cause infection, product failure, overheating leading to burn injuries. (Refer to "9. Post-use Maintenance")
- Do not change or otherwise reinsert a bur into the handpiece during a procedure. Doing so can introduce blood or other debris into the chuck, increasing the risk of cross-contamination.
- Do not rotate the motor when the Bur Lock Ring is at the OPEN position, or when a bur is not mounted in the chuck. Either situation It may cause burn injuries due to sudden heat generation. In addition, the motor insert could be damaged which may lead to difficulty in connecting and disconnecting the handpiece.
- Always keep the bur shank clean. Dirt or debris in the chuck could cause heating, poor bur concentricity, or low chuck retention force.
- Supply coolant water through the irrigation nozzle while using the handpiece. Other irrigation methods than recommended may lead to overheating causing burn injuries.
- Do not apply excess pressure to the bur as it may break, bend or become difficult to remove. Doing so may also cause burn injury due to abnormal heating.
- Ball bearings are wear items and may require replacement. Be sure to inspect the handpiece before use (Refer to "8. Checking the Handpiece Before Each Use"). If abnormal vibration or noise are found, the bearings or other internal parts may be worn out or damaged. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- The temperature at the tip of the handpiece may reach up to 51°C depending on the environment.
Contact with the body surface of the patient or the operator must be within 1 min..
Contact for longer than the specified time could cause higher risk of low temperature burn injury.

CAUTION

- When operating the product always consider patient safety.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer. (Refer to "8. Checking the Handpiece Before Each Use")
- Do not disassemble or alter the handpiece except as recommended by NSK in this Operation Manual.
- Do not allow any impact to the handpiece. Do not drop the handpiece. Deformation may cause the handpiece to fail during use.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur.
 - Bur which is cracked on the edge or axis.
 - Tampered bur.
- Do not use burs other than specified in "11. Specifications" as such burs may be accidentally released or break during rotation.
- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the handpiece head bearings, or accidental release of the bur, may occur.
- Always follow the instructions provided by the bur manufacturer.
- Always confirm before use that the bur is attached securely. Improper attachment may lead to sudden release of the bur during rotation.
- Never turn the Bur Lock Ring while the handpiece is rotating.

CAUTION

- For patients safety, be sure to use a sterile handpiece and sterile bur for each surgical operation and / or dental procedure.
- Should it be necessary to change or replace a bur during treatment, make sure that clean and sterile handpiece and bur are used.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Do not use the following fluids to wipe, immerse or clean the product; strong / super acid water, strong acid / alkaline chemicals, chlorinecontaining solutions, solvents such as benzine or thinner.
- The handpiece is delivered in a non-sterile condition and perform steam sterilization prior to use.
- Perform regular function and maintenance checks. (Refer to "10. Periodical Maintenance Checks")
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

NOTICE

- Contact your Authorized NSK Dealer if repairs are necessary.

3. Accessory List

No.	Part Name	Quantity
1	Irrigation Nozzle	1
2	E-Type Spray Nozzle	1
3	Test bur	1

4. Connecting and Disconnecting the Handpiece from the Motor

4-1 Connecting

- 1) Insert the handpiece until it reaches the end of the Motor Insert. (Fig. 1)
- 2) For motors that feature a locking mechanism, insert the handpiece until it reaches the end of the motor insert and then screw the handpiece until it is locked into place.
- 3) After connection, confirm that the handpiece is firmly connected to the motor by pulling and pushing the handpiece.


Fig. 1

4-2 Disconnecting

Hold the Motor and the Handpiece, then pull apart.

CAUTION

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964) with the Max. Speed of 40,000 min⁻¹ or less.

5. Mounting and Removing the Bur

5-1 Mounting the Bur

- 1) Turn the Bur Lock Ring in the direction of the "🔓" mark until it clicks. (Fig. 2)
- 2) Insert the bur fully into the chuck then turn the Bur Lock Ring in the direction of the "🔒" mark until the "●" mark and the "🔒" mark align with a click.
- 3) Confirm that the bur is secure by gently pulling and pushing the bur.


Fig. 2

5-2 Removing the Bur

Turn the Bur Lock Ring in the direction of the "🔓" mark until it clicks then remove the bur.

WARNING

- Always use a clean, sterile bur. Be sure to wear sterile gloves before mounting the bur.
 - Do not rotate the motor when the Bur Lock Ring is at the OPEN position (the "●" mark and the "🔓" mark are misaligned), or when a bur is not mounted in the chuck. Either situation It may cause burn injuries due to sudden heat generation. In addition, the motor insert could be damaged which may lead to difficulty in connecting and disconnecting the handpiece.
 - After mounting the bur in the chuck, be sure to confirm that the "●" mark and the "🔓" mark align. If these marks are not completely aligned, burn injuries may occur due to handpiece overheating during rotation.
 - Do not replace, remove / mount the bur during treatment. This causes blood or debris to enter the handpiece, leading to secondary infection. If you need to use a different bur during treatment, mount it on a new handpiece.
 - Do not change or otherwise reinsert a bur into the handpiece during a procedure. Doing so can introduce blood or other debris into the chuck, increasing the risk of cross-contamination.
- Should it be necessary to change or replace a bur during treatment, make sure that clean and sterile handpiece and bur are used.

CAUTION

- Do not mount or remove the bur until the motor has completely stopped.
- Insert a sterilized bur and always keep the bur shank clean. Should it be necessary to change or replace a bur during treatment, make sure that clean and sterile handpiece and bur are used.

6. The Bur

The SGS / SGA series handpieces are designed for use with surgical burs (Ø2.35mm, 70mm in shank length).

■ To Mount the Stopper (ONLY for SGS-ES / SGA-ES)

When using a bur designed for dental straight handpieces in implant surgery, etc., attach the Stopper (Optional) to the handpiece. The Stopper is designed for burs with the following specifications:

Ø2.35mm, 44.5mm in total bur length, 30mm or more in chuck length. (Fig. 3)

- 1) Remove the bur from the handpiece. (Refer to "5. Mounting and Removing the Bur")
- 2) Turn the Bur Lock Ring to the OPEN position. (Fig. 2)
- 3) Insert the threaded portion of the Stopper into the chuck then screw it in with the provided screwdriver.

*Refer to Option Parts List to identify the correct parts.

CAUTION

- Use only with surgical burs designed for the following dental straight handpieces.
SGS-E2S, SGA-E2S


Fig. 3

7. Irrigation Nozzle

Insert the Irrigation Nozzle (Accessory) into the handpiece. Then push the Irrigation Tube into the nozzle securely. (Fig. 4)


Fig. 4

8. Checking the Handpiece Before Each Use

Wear sterile gloves then follow the check below before use. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 1) Mount the Bur. (Refer to "5. Mounting and Removing the Bur")
- 2) Adjust the handpiece rotation speed in accordance with the allowable maximum speed for the attached bur.
At that time, do not exceed the Max. rotation speed for the handpiece. (Refer to "11. Specifications")
- 3) Rotate the handpiece and check for abnormal rotation, vibration, noise and overheating.
Also make sure that the water spray is hitting the tip of the bur.
If any of these conditions are found, stop using the handpiece and contact your Authorized NSK Dealer.
- 4) Check that water is supplied to the cutting portion of the bur. If any abnormalities are found with the Irrigation function, stop using the Irrigation Nozzle and replace it with a sterilized new Irrigation Nozzle. (Refer to "15. Spare Parts List")
- 5) If there are no immediate apparent abnormalities, continue to rotate the handpiece for about 1 minute. Then stop rotation and check for abnormal heating at the tip of the handpiece. If abnormal heating is detected, stop using the handpiece and contact your Authorized NSK Dealer.
- 6) Operate the handpiece if no abnormalities are observed.

CAUTION

- To avoid injury, keep your hands away from the bur during rotation.

9. Post-use Maintenance

Post-use maintenance includes cleaning, lubricating, and sterilization after each patient use.

WARNING

- To ensure cleaning and sterilization efficacy, use only the following procedures for reprocessing.
- Follow any additional local directives, standards, and guidelines for cleaning, disinfection, and sterilization.
- Clean, disinfect, lubricate and sterilize the handpiece immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the handpiece may cause infection, product failure, overheating leading to burn injuries.

9-1 Cleaning at point-of use

CAUTION

- Do not use the following liquids to wipe, immerse, or clean the product:
strong / super acid water, strong acid / alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- Do not immerse NSK instruments in disinfectant solution or clean in ultrasonic devices.

- 1) Always wear protective gloves, a mask, and protective goggles for safety purposes and to minimize the risk of infection. (Fig. 5)


Fig. 5

- 2) Operate the handpiece while running clean water into the irrigation tube with the irrigation nozzle, at the chair-side for at least 20 seconds to purge fluids out of the handpiece. (Fig. 6)
After the handpiece rotation has completely stopped, confirm the handpiece head is NOT heating abnormally. (Water temperature 38°C or less, drinking water level of water quality recommended)


Fig. 6

⚠ CAUTION

- If the head is heating abnormally, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Be sure to take the step "9-1 2)" before removing the bur. If the handpiece operates without the bur, it may cause burn injuries due to sudden heat generation. In addition, the motor insert could be damaged which may lead to difficulty in connecting and disconnecting the handpiece.

3) Remove the bur from the handpiece using tweezers to firmly hold the bur. (Fig. 7)


Fig. 7

- 4) Wipe the exterior of the handpiece clean using a cloth moistened with disinfectant or ethanol with a concentration of 60 - 90%. (Fig. 8)

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada.

CaviCide[®], CaviWipes[®] (manufactured by Metrex)


Fig. 8

- 5) Remove the irrigation tube from the irrigation nozzle and remove moisture with a dry cloth. Then remove the irrigation nozzle from the handpiece. (Fig. 9)


Fig. 9

6) Remove the handpiece from the motor. (Fig. 10)
(Refer to "4-2 Disconnecting")

7) Carry the handpiece and the irrigation nozzle to the decontamination area.


Fig. 10

9-2 Cleaning, Disinfecting, and Drying

Perform cleaning, disinfecting, and drying using either the manual or automatic method.

Manual Method (Cleaning, Disinfecting, and Drying the Exterior)

Wash the handpiece and irrigation nozzle under running water following the procedures below.
(Water condition: $\leq 38^{\circ}\text{C}$, $\geq 3.5\text{L} / \text{min}$; water should be of the same quality as drinking water)

- 1) Clean the external surface of the handpiece with a soft bristled toothbrush for 15 seconds or more.
(Fig. 11)


Fig. 11

- 2) Clean the interior cavity of the handpiece for 15 seconds or more, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409. (Fig. 12)


Fig. 12

- 3) Clean the external surfaces of the irrigation nozzle with a soft bristled toothbrush for 15 seconds or more. (Fig. 13)


Fig. 13

- 4) Clean the attachment part of the irrigation nozzle for 15 seconds or more, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409. (Fig. 14)


Fig. 14

Under appropriate lighting (natural office lighting, 500lux or higher), inspect the handpiece and irrigation nozzle for blood or other debris. If any visible debris remains, repeat the process until the handpiece is visually clean.

5) Wipe the exterior of the handpiece and the irrigation nozzle clean using a cloth moistened with disinfectant or ethanol with a concentration of 60 - 90%. (Fig. 15, 16)

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada.

CaviCide[®], CaviWipes[®] (manufactured by Metrex)

After cleaning and disinfection, if the residual moisture is present, wipe off with dry cloth or blow it off with compressed air($\leq 0.35\text{MPa}$) until there is no moisture in the interior and exterior.

When blowing off with compressed air, cover the handpiece with cloth to prevent scattering of water.

Proceed to "9-3 Cleaning and Lubrication".


Fig. 15


Fig. 16

Automatic Method (Cleaning, Disinfecting, and Drying the Exterior and Interior)

NSK handpieces with this symbol are compatible with medical thermo-disinfectors.

Due to the variation in cleaning / disinfection cycles and cleaning agents available from the various manufacturers, NSK has validated the Getinge 46 Series Washer Disinfector Cycle P7-Anesthesia for automatic processing / reprocessing of the product.

Only validated procedures must be used for cleaning, disinfection and sterilization. When using reprocessing procedures that are different from those described in this manual, those procedures must be validated by the respective practice or hospital using a thermo-disinfector that conforms to ISO 15883-1.

NOTICE

- All water is "building supply" (at least drinking water quality) unless otherwise indicated.

- 1) Prepare and pre-disinfect the device as described above in step "9-1 Cleaning at point-of use".
- 2) Wipe off a dirt on a part held by a handpiece holder using the cloth moistened with 60 - 90% ethanol or disinfectant.
When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada.

CaviCide[®], CaviWipes[®] (manufactured by Metrex)

3) Perform main wash under the following conditions. Use a Getinge 46 Series thermo-disinfector (or equivalent device):

Pre-Wash 1: Cold water (<110°F (43°C)) for 2 min.

Pre-Wash 2: None - Not Applicable

Enzyme Wash: Hot tap water and 60ml of heated (95°F (35°C)) detergent (Getinge Clean Renuzyme) for 3 min. at 120°F (49°C)

Detergent Wash: None - Not Applicable

Rinse 1: Hot tap water for 2 min.

Rinse 2: Hot tap water for 2min.

RO Final Rinse: Hot tap water for 32 min. at 167°F (75°C)

Drying: 20 min. at a temperature of 140 - 194°F (60 - 90°C)

Upon removal from the thermo-disinfector, if the residual moisture is present, wipe off with dry cloth or blow it off with compressed air ($\leq 0.35\text{MPa}$) until there is no moisture in the interior and exterior.

4) After cleaning, disinfecting, and drying the handpiece, confirm that it is clean under appropriate lighting (500lux or higher). If dirt still remains, repeat the procedures.

5) Proceed to "9-3 Cleaning and Lubrication".

CAUTION

- Handpieces must be removed from the thermo-disinfector immediately (within 1 hour) after the cleaning, disinfecting and drying cycle is complete to prevent corrosion.
- After washing with thermo-disinfector, prior to lubrication, dry the product until all internal moisture is completely removed. Thermo-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

9-3 Cleaning and Lubrication

Perform cleaning and lubrication using the manual method.

Manual Method (Cleaning and Lubricating the Interior)

- 1) Check that the " ● " mark and the " 🔒 " mark align.
- 2) Attach an E-Type Spray Nozzle to the nozzle of the PANA SPRAY Plus.
- 3) Shake the can 3 or 4 times vigorously, and insert the E-Type spray nozzle into the rear of the handpiece while holding the head of the handpiece with a cloth etc..
- 4) Hold the handpiece, and spray for 2 or 3 seconds until oil comes out of the handpiece head (Fig. 17).

Repeat lubrication until foreign material stops coming out of the tip.


Fig. 17

⚠ CAUTION

- Check that the " ● " and the " 🔒 " mark of the Bur Lock Ring are aligned. Spraying with the " ● " and the " 🔒 " aligned will lead to increased pressure of the spray, which could cause the handpiece to eject from the spray.
- Do not use sprays other than NSK PANA SPRAY Plus. Using sprays other than PANA SPRAY Plus may cause the handpiece to overheat.
- Hold the spray can upright.
- Firmly hold the handpiece to prevent it slipping when spray pressure is applied.
- Spray lubricant until it expels from the handpiece head.

NOTICE

- NSK recommends the use of "Spray Mist Absorber" (Y900084) to prevent oil mist expelling out of the handpiece head.

9-4 Packaging, Sterilizing, and Drying

- 1) Insert the handpiece into an FDA-cleared sterilization pouch that conforms to ISO 11607-1 and seal the pouch.
- 2) After packaging in a sterile pouch, sterilize using the following validated cycles. Use an FDA-cleared steam sterilizer to perform sterilization. Follow local rules, regulations, and guidelines regarding the reprocessing of devices.

Type	Gravity Displacement	Pre-Vacuum (Dynamic Air Removal)
Temperature	132°C	132°C
Full Cycle Time	15 min.	4 min.
Drying Time	30 min.	30 min.

CAUTION

- Use an FDA-cleared steam sterilizer to perform sterilization.
- Do not touch the product immediately after steam sterilization as it will be very hot and must remain in a sterile condition.
- Do not perform steam sterilization on the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.

CAUTION

- Clean and lubricate the handpiece prior to sterilization. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- Be sure to use sterilizers that can perform sterilization up to 135°C. In some sterilizers, the chamber temperature may exceed 135°C. Do not use these sterilizers as failure of the handpiece could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.

9-5 Storage

Store the product in a dry, clean location.

CAUTION

- After the sterilization and drying cycles are complete, remove the handpiece immediately from the sterilizer to store it.
- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in "11. Specifications".
- Sterilization is not guaranteed after the sterilization retention period specified by the manufacturer and seller of the sterilization pouch has elapsed. If the sterilization retention period has elapsed, perform sterilization again with a new sterilization pouch.

10. Periodical Maintenance Checks

Every 3 months perform periodical maintenance checks, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Bur mounting / removing	Make sure that the bur can be properly mounted / removed. After mounting, make sure that the bur is firmly fixed by pulling and pushing it.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.

11. Specifications

Model	SGS-ES	SGS-E2S	SGA-ES	SGA-E2S
Handpiece Type	Straight		20° Angle	
Allowable Max. Speed	40,000min ⁻¹	80,000min ⁻¹	40,000min ⁻¹	80,000min ⁻¹
Max. Rotation Speed (Motor)	40,000min ⁻¹			
Max. Rotation Speed (Handpiece)	40,000min ⁻¹	80,000min ⁻¹	40,000min ⁻¹	80,000min ⁻¹
Gear Ratio	1 : 1 Direct Drive	1 : 2 Increasing	1 : 1 Direct Drive	1 : 2 Increasing
Bur Type	Short Surgical Burs Ø2.35 x L70 mm			
Chuck Length	55.5mm	53mm	55.5mm	53mm
Max. Working Part Diameter	Ø7.0mm	Ø3.0mm	Ø7.0mm	Ø3.0mm
Water Pressure	Min. 50ml / min (0.2MPa (2.0bar))			
Use Environment	Temperature: 10 - 35°C, Humidity: 30 - 75% (No Condensation)			
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Pressure: 500 - 1,060hPa			

12. Symbol

Open

Close

Position

This product can be sterilized in a steam sterilizer up to Max. 135°C.

This product can be washed via Thermo Disinfector.

Conforms to CE European Directive of "Medical device directive 93 / 42 / EEC."

Manufacturer.

Consult operating instructions.

Caution: U.S.Federal law restricts this device to sale by or on the order of a licensed physician.

13. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

14. Option Parts List

Model	Order Code	Compatible Product
Stopper	Y900160	SGS-ES / SGA-ES

15. Spare Parts List

Model	Order Code
E-Type Spray Nozzle	Z019090
Irrigation Nozzle	H263018

16. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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NSK America Corp www.nskdental.com